

Cryotainer B.V.
Fokkerstraat 507
3125 BD Schiedam
Phone: +31 (0)10 262 30 53

P.O. Box 130
3100 AC Schiedam
www.cryotainer.com



Datasheet Twin Ambient Vaporizer Unit



| | | |
|-------------------------------------|--------------------|--|
| Dimensions l x h x b | mm | 12.200 x 3.650 x 2.500 |
| Maximum working pressure | bar | 40 |
| Maximum flow | Nm ³ /h | 3.400 (4 x 850) for 8 hours per day maximum |
| Maximum outlet temperature nitrogen | °C | 15 ⁰ till 20 ⁰ below ambient temperature |
| Total weight | kg | 38.000 |
| Suitable for cryogenic gasses | | LIN; LOX*; LAR; LCO ₂ |
| Outlet connection | | Four sides 2" weco fig. 100 |

* For use of LOX container and vaporizer should be purged before use